

INSTRUCTION MANUAL

FOOTSTOOL



First of all we wish to thank you for the trust you placed in us by selecting a **VERMEIREN** product.

Before using the product, read the instruction manual carefully: it will familiarise you with this product.

You will also be given hints about the care and service of your footstool. Please remember that observance of our hints would keep your product in an excellent condition of preservation and it will function perfectly even after years of use.

If you have any further questions, please consult your specialist dealer. **To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility.**

TECHNICAL DETAILS

Height:	460 mm
Width:	490 mm
Length:	480 mm
Minimum angle:	-25°
Maximum angle:	0
Weight:	5 kg

We reserve the right to introduce technical changes. Measurement tolerance +/- 15 mm / 1,5 kg

CONTENTS OF CONSIGNMENT

Unpack your product and verify if the delivery is complete. Following items should be included:

- Footstool
- Instruction manual

Verify your product for transport damage. If you find any damages after delivery proceed as follows:

- contact the transporter
- have a list made of any problems
- contact your supplier

ASSEMBLY INSTRUCTIONS

The footstool is delivered fully assembled. Place the footstool on level surfaces only and make sure that all four legs stand solidly on the floor.

ADJUSTING ANGLE

Pull the footstool on the front side upwards to the desired height.



To lower the surface, pull the footstool upholstered surface to the upper position. The rastomat will then remove the lock. You can now completely lower the upholstered surface.

SAFETY

- ➔ When transporting the footstool, carefully guard against bruising and shearing - danger of injury!
- ➔ Before using, check whether all four legs stand evenly on the floor – risk of tipping over the footstool!
- ➔ Repairs and disinfection should only be done by authorised persons !
- ➔ Do not place your fingers between the angle adjustment system.
- ➔ Do not stand on the footstool.
- ➔ Risk of burns – Be careful when using in hot or cold environments (sunshine, extreme cold, saunas, etc.) for a sufficient amount of time and when touching.
- ➔ Take note of the instructions for care and service. The manufacturer is not liable for damage caused by improper servicing/care.

CARE

Cleaning shall be according following instructions:

- **Upholstery / frame**
 - Clean upholstery / frame with a cloth moistened with hot water. Be aware that you do not soak the upholstery / frame.
 - Use a mild commercial detergent for removing stubborn dirt.
 - Stains can be removed by using a sponge or a soft brush.
 - Do not use strong cleaning liquids like solvents, nor use hard brushes.

- Never clean with steam and/or pressure washers.
- The manufacturer is not liable for damage caused by improper cleaning agents.

- **Surface layer**

The high quality of the surface layer guarantees optimal protection against corrosion. If the outer coat is damaged by scratching or in some other way, get your specialist dealer to repair the affected surface.

When cleaning, only use warm water and normal household detergents and soft brushes and cloths. Ensure that no wetness gets into the tubes.

Initially the zinc parts only require rubbing with a dry cloth. Stubborn dirt can best be removed by using a suitable commercial zinc polish.

SERVICE

Lifetime of the footstool is influenced by its use, storage, regular maintenance, servicing and cleaning.

The footstool must be checked at least every second year by authorised persons/specialist dealers, and at least every time before resuming its use. The following parameters must be checked:

- Check the frame for plastic deformation and reduced functioning.
- Visual inspection of damage to the paintwork (danger of corrosion)
- Check the condition of the cushions (tears, brittleness, soaking)
- Completeness of consignment ?
- Instruction manual available?
- Cleaning

The service should only be confirmed in the "Service plan" if at least the above aspects have been checked. If your specialist dealer does not undertake the services which you require, kindly contact the manufacturer. We will gladly refer you to authorised specialist dealers in your vicinity.

DISINFECTION

All parts of the footstool can be disinfected by scrubbing. In principle all surfaces of a system or a product have to be disinfected before handing it over to another user (re-use).

➔ Only authorised technicians trained in the functioning and application of disinfectants should apply disinfectants.

➔ Wear suitable protective clothing to prevent skin irritations which might be caused by the disinfectant. Also note the product information of the solutions used.

➔ The manufacturer does not accept liability for damage and injury caused by improper implementation of the disinfection.

Use only disinfectants included in the list of the Robert Koch Institute and the verified procedures of disinfection by scrubbing or wiping.

The current situation in regard to the disinfectants included in the RKI list, can be obtained from the Robert Koch Institute (RKI, home page: www.rki.de).

SHIPPING AND STORAGE

The shipping and storage of the footstool shall be according following instructions:

- Store in a dry place (between +5°C and +41°C).
- The relative humidity of the air should be between 30% and 70%.
- Provide sufficient covering or packaging to protect the toilet chair from rust and foreign bodies (e.g. salt water, sea air, sand, dust).

WARRANTY

Excerpt from the "General Business Conditions".

(...)

5. The guarantee period for warranty claims is 24 months.

(...)

The warranty excludes damage arising from structural changes to our products, insufficient maintenance, defective or improper handling or storage or the use of non-original parts.

DISPOSAL

When disposing of the footstool, contact your local disposal centre or return the product to your specialist dealer who, after submitting it to a hygienic procedure, will be able to send it back to the manufacturer who will dispose of and recycle it correctly, separating it into its component materials.

Packaging materials can be taken to disposal or recycling centers or to your specialist dealer.

SERVICE PLAN

The footstool was serviced:

Dealer's stamp

Date :

Dealer's stamp

Date :

Dealer's stamp

Date :

Dealer's stamp

Date :

For service checklists and additional technical information, please see our specialist dealers nearest to you. More information on our website at:

www.vermeiren.com

www.vermeiren.be



VERMEIREN GROUP

Vermeirenplein 1/15
2920 Kalmthout

Tel.: +32 (0)3 620 20 20

Fax: +32 (0)3 666 48 94

Website: www.vermeiren.be

E-mail: info@vermeiren.be