


Bathing seats

EASY LIFE RANGE



All our *Easy life* products are ideal for usage in the shower, bathtub, bathroom or other humid environments.

The special designed rubbers at the end of each leg will create a safe stability for the bench and avoid the sliding of the bench.



Jill



Grace

Removable shower seat for the wall.
Easy installation.



Jinny



Lory

Bathroom chair

Both Lola and Lory are two practical, height adjustable bath stools with rubber feet.

All 4 legs can be adjusted in height by simply pressing a button. For easy storage, all 4 legs can be slid inwards.



Sidney

Bathroom seat with toilet opening



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

