

Pratic

Electrical adjustable homecare bed
EN 1970 / A1



Movable lifting pole,
equipped with a
double serumholder.
Can be placed left or right.

**AVAILABLE WITH
2 or 3 FUNCTIONS
SIDESUPPORTS in WOOD**










Linak IP66 motor works ste-
pless and silent. Easy to
disassemble without need
for extra tools.

Standard equipped
with four wheels with
separate brakes.



Technical information Pratic bed

 width	 length	 min. position	 max. position	 inclination head part	 inclination foot part	 users weight
1030	2080	360	790	83°	30°	140 kg

Sizes in mm, weight in kg

Deviations in sizes and weight (1,5 mm and 1,5 kg) are possible

Comfortable and practical

The Pratic bed comes with an electrical adjustable headsupport. One push on the button will place you in the ideal position to read, watch tv or dine. The foot-part can be placed in your preferred position, either manual or electrical (option).

This high-low bed comes with four wheels with separate brakes which enhances its stability. The both strong and flexible slats are easy to maintain. Furthermore, the bed is easy to assemble without any need for extra tools.

Esthetic and maintenance friendly

The wooden front and back panels provide this high-low bed with a homely character. The warm wooden appearance changes this bed in a fine looking piece of furniture that fits into any environment. On top of that, the bed is very easy to maintain.



Colors

Frame	Panels
	
White	Beech



Comes standard with a transporterset

Options

- Centrale brake



Certificate TÜV SÜD
EN 60601-1/A2:1995
EN 60601-2-38/A1:2000
EN 1970/A1:2005

N.V. Vermeiren N.V.

Vermeirenplein 1-15
 B-2920 Kalmthout
 Belgium

Tel.: +32 (0) 3 620 20 20
 Fax: +32 (0) 3 666 48 94

Visit our website: www.vermeiren.com

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility