

Mercurius 4D

Stylish and timeless

A stylish and fully functional 4-wheel scooter. Sporty rims, modern shopping basket, two rear view mirrors, fully adjustable and swivel seat with comfortable headrest and adjustable armrests.

LED lighting at the front and rear as well as direction indicators ensure maximum safety even in poor visibility. Also equipped with two USB ports to charge your devices on the go.



Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility

Technical data Mercurius 4D

Total width	680 mm
Total length	1430 mm
Total height	1355 mm
Seat height	470- 495 mm
Turning radius	3740 mm
Batteries	2x12V - 80Ah

Driving range	ca. 40 km
Motor	1000 W
Speed	15 km/u
Weight	116,75 kg
Users weight	150 kg
Additional information	CE

Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Theoretical range, calculated under ideal conditions in accordance with EU standards. Speeds and range may deviate (approx. 15%).



Options

- Rollator holder
- Cover
- Luggagebox with lock
- Crutch holder



Colors



Your guarantee of quality

Regarding quality, the **Mercurius 4D** complies with the highest norms, laid down by certified test institutes in most European countries.

N.V. VERMEIREN N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility

Subject to technical alterations