

Mercurius 4 Limited Edition

- 4 wheels
- Length: 1430 mm
- Width: 680 mm
- Height: 1355 mm
- Action radius: ca. 45 km*
- Speed: 15 km/h
- Batteries: 2x12 V / 75 Ah
option 2x12 V / 80 Ah
- Users weight: 150 kg
- Color: Champagne

outdoor scooter



We reserve the right to introduce technical changes. Measurement tolerance ± 15 mm / 1,5 kg / °.
* Theoretical driving distance, calculated under best conditions in accordance with the European Standards.

VERMEIREN

Opening up new horizons in mobility

Vermeirenplein 1-15
B-2920 Kalmthout - Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

info@vermeiren.be
www.vermeiren.com



N.V. Vermeiren N.V.