

Quava

ROLLATOR

total width	655 mm
seat depth	225 mm
seat width	450 mm
seat height	540 mm
total height	830 - 955 mm
total length	710 mm
weight with bag	7,05 kg
max. users weight	136 kg
wheels (front/rear)	200 x 35 mm
colours	red - carbon grey

Measurement tolerances: ± 1.5 mm / 1,5kg.
We reserve the right to introduce technical changes.



Model shown only for illustration purposes. Certain features may be optional.

N.V. Vermeiren N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: +32 (0)3 620 20 20 • Fax: +32 (0)3 666 48 94 • info@vermeiren.be • website: www.vermeiren.com



Opening up new horizons in mobility