

# V500 30°

Reclinable comfort

The V500 30° is a light weight manually propelled wheelchair for temporary or permanent use with maximized settings. Similar to the V500 but with a back inclination up to 30°. Complete with special designed armrests which are adjustable, foldable and removable.



**VERMEIREN**  
we care for you

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.



Push bar



Anti-tipping wheels



Tension adjustable backrest



Detail brake

## Quality guaranteed

Regarding quality, the **V500 30°** complies with the highest norms, laid down by certified test institutes in most European countries.

## Colours



## Technical data

■ Total width	590/620/640/660/680/700 mm
■ Total height	860 - 920 mm
■ Total length	975 - 1035 mm
■ Folded width	310 mm
■ Seat width	390/420/440/460/480/500 mm
■ Height armrests	150 - 210 mm

■ Seat height	440 - 530 mm
■ Height backrest	440 mm
■ Seat depth	380 - 500 mm
■ Weight	16,6 kg
■ Maximum user weight	130 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you