

# 708D

## Standard manual wheelchair

**Aesthetics, reliability and functionality. All removable parts such as footrests, armrests, wheels and cushions can be replaced thanks to the modular system.** The wide range of accessories and options allows the wheelchair to be optimally adapted to the user. Well known and recognized for its sturdiness and quality, the 708D makes an ideal asset for a rental fleet or as a durable choice to own for a very long time.



**VERMEIREN**  
we care for you

## Features



**Aluminium footplates complete with rubber protection wheels**



**Easy to fold, easy to transport**



**Hemiplegic (HEM2) version. Double rims to propel and to steer**



**Amputee (AMP) version. Position of the rear wheel has shifted to the back**

- 3 types of folding and removable armrests with long and comfortable armrests
- Wide range of options and accessories
- Armrest locking system is easy to access and to operate
- Various seat and backrest options
- Available in 6 seat widths
- 4 types of footrests and 5 types of leg rests


## Colours



C29  
grey

## Technical data

Total width	560 / 590 / 610 / 630 / 650 / 670 mm
Total height	910 mm
Total length	1109 mm
Folded width	280 mm
Seat width	390 / 420 / 440 / 460 / 480 / 500 mm
Height armrests	225 mm

Seat height	490 mm
Height backrest	420 mm
Seat depth	430 mm
Weight	18 kg
Maximum user weight	120 kg
Additional info	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **708D** complies with the highest norms, laid down by certified test institutes in most European countries.

## VERMEIREN GROUP

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you