

253

ROLLATOR

total width	560 mm
seat width	415 mm
seat height	530 mm
total height	795 - 915 mm
total length	670 mm
max. users weight	130 kg
total weight	4,60 kg
foldable	yes
folded width	185 mm
removable basket	yes
removable plexi tablet	yes
frame color	inox

Deviations in sizes and weight (15 mm and 1,5 kg) are possible - Subject to technical alterations.



N.V. Vermeiren N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: 00 32 (0)3 620 20 20 - Fax: 00 32 (0)3 666 48 94 - info@vermeiren.be - visit our website: www.vermeiren.com