

692

Foldable ramps

Placing a scooter or a wheelchair in a car (for transport) becomes easy thanks to these ramps. These foldable ramps make sure that you can safely place heavy scooters and wheelchairs in the car.




**Including storage bags.
One set contains two ramps.**



Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility

Technical data 692



Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible

01 Handles to easily transport the ramps

02 Non-slip surface for better traction of the wheelchair or scooter

Maximum load **03**
150 kg per ramp

Using the ramps

The ramps are easy to use by performing the following steps:

1. Unfold the plate.
2. Place the anti-slip overlap on the surface to be covered.
3. Ensure the ramps are stable.
4. Without occupant (!), ride the scooter or wheelchair on the ramps.

04 Provided with non-slip overlay for better grip of the ramps



N.V. VERMEIREN N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility



Subject to technical alterations