



D200-V

Light and practical

The D200-V is a light and multi-functional wheelchair.

A modern and striking blue frame color combined with several adjustments to personalize this wheelchair to your own configuration. Foldable in mere seconds (if needed the armrests, legrests and wheels can be detached) and your D200-V is ready for transport.



VERMEIREN
we care for you

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **D200-V** complies with the highest norms, laid down by certified test institutes in most European countries.

Colours

C80

blue



Quick-Release rear wheel



Removable legrests



**Option B78
Anti-tipping wheels**



**Option BZ8 - Comfort and
length adjustable legrests**



**Push handles
adjustable in height**



**Armrests adjustable
in height and depth**



**Tension adjustable backrest
Height adjustable backrest**



**Adjustable inclination
of the footplates**

Technical data

■ Total width	580 / 600 / 620 / 640 / 660 / 680 / 700 / 720 mm
■ Total height	890 - 1060 mm
■ Total length	995 - 1045 mm
■ Folded width	300 mm
■ Seat width	380 / 400 / 420 / 440 / 460 / 480 / 500 / 520 mm
■ Height armrests	200 - 240 mm

■ Seat height	470 - 530 mm
■ Height backrest	450 - 510 mm
■ Seat depth	395 - 450 mm
■ Weight	15,8 kg
■ Maximum user weight	130 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer