



# V500 Light

**OUR LIGHTEST  
STANDARD WHEELCHAIR**



**VERMEIREN**  
we care for you



## Our lightest

Lighter than the V500, the V500 Light is our lightest standard wheelchair. Whether for temporary or permanent use, it improves the user's daily life and makes transport and transfers easier.

## Highly configurable

As part of the V-series, the V500 Light offers a wide range of accessories and options to create and modify your V500 Light.



*Toolfree removable armrests*

*In and outward folding footrests  
Removable footrests*

*Folding brakes levers for easy transfers*

*Foldable for easy transport or storage*



*Adjustable seat height (4 positions)*

*Adjustable seat depth: 42 to 48 cm*

*Height adjustable push handles and backrest*

*Adjustable seat angle: 0° - 10°*

*Adjustable armpads in height and in depth*



# V500 Light

## Available colors



C80



C86

## More info



## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

The **V500 Light** offers a wide range of accessories and options. Consult our online orderform to get a complete overview.

## Technical data

■ Total width	590 - 620 - 640 - 660 - 680 - 700 mm
■ Total height	860 - 920 mm
■ Total length	975 - 1035 mm
■ Folded width	310 mm
■ Seat width	390 - 420 - 440 - 460 - 480 - 500 mm
■ Height armrests	150 - 210 mm

■ Seat height	470 - 530 mm
■ Height backrest	400 - 420 - 440 - 460 - 480 mm
■ Seat depth	420 - 480 mm
■ Weight	15,5 kg
■ Maximum user weight	130 kg
■ Additional info	<b>CE</b>

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group N.V.

Vermeireplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you