



# D200-V

Light and practical

**The D200-V is a light and multi-functional wheelchair.**

A modern and striking blue frame color combined with several adjustments to personalize this wheelchair to your own configuration. Foldable in mere seconds (if needed the armrests, legrests and wheels can be detached) and your D200-V is ready for transport.



**VERMEIREN**  
we care for you

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **D200-V** complies with the highest norms, laid down by certified test institutes in most European countries.

## Colours

**C80**

blue



**Quick-Release rear wheel**



**Removable legrests**



**Option B78  
Anti-tipping wheels**



**Option BZ8 - Comfort and  
length adjustable legrests**



**Push handles  
adjustable in height**



**Armrests adjustable  
in height and depth**




**Tension adjustable backrest  
Height adjustable backrest**



**Adjustable inclination  
of the footplates**

## Technical data

Total width	580 / 600 / 620 / 640 / 660 / 680 / 700 / 720 mm
Total height	890 - 1060 mm
Total length	995 - 1045 mm
Folded width	300 mm
Seat width	380 / 400 / 420 / 440 / 460 / 480 / 500 / 520 mm
Height armrests	200 - 240 mm

Seat height	470 - 530 mm
Height backrest	450 - 510 mm
Seat depth	395 - 450 mm
Weight	16,2 kg
Maximum user weight	130 kg
Additional info	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN GROUP

Vermeirenpolein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer