

Never leave without your canopy



Protected in all weather conditions

The ultimate 4 seasons canopy for mobility scooters



VERMEIREN
we care for you

Canopy for mobility scooters. REDESIGNED.



No matter the weather, the **Movo canopy** protects you first of all from the sun, but also from all weather conditions. Combining protective fabric against the rain, windshield and sun visor, it also protects against cold and wind.

Easy to assemble, easy to install, and of impeccable design quality, the Movo canopy protects you from the sun, rain, cold and wind, while maintaining optimal ventilation in the cabin

Maximize the use of your vehicle

A much longer season

Get your vehicle out earlier in the spring, and extend your mobility scooter season, thanks to its protection against cold and wind.

Double your annual outings Increase your number of annual outings in uncertain weather. Don't be afraid of getting caught in the rain far away from home thanks to the retractable canopy, which can be adapted according to the weather.

A canopy that turns into a cabin in an instant

Thanks to its retractable technology, it does not need to be lowered at all times. You can choose the position that suits you, depending on the weather. You therefore have the option: canopy and cabin option with the same product.



Super-fast installation

You can install or remove the Movo canopy in less than 10 seconds. It's ultra-quick.



Easy to assemble

A detailed assembly manual and an accurate assembly video are available on our website, so you will be able to assemble the canopy in less than 15 minutes. Once assembled, the installation can be completed in less than 10 seconds.



Compact design

The width of the Movo canopy is 74 cm, which allows it to fit through all openings, even small interior doors.



Ultralight

With a total weight of 8,7 kg (19 lbs), the Movo canopy can be installed easily by one person.

No matter the weather

The Movo canopy protects you first of all from the sun, but also from all weather conditions. Combining protective fabric against the rain, windshield and sun visor, it also protects against cold and wind.



For your scooter

Thanks to an ingenious installation using the hitch under the seat, Movo canopy is compatible with most models of **VERMEIREN** mobility scooters*

* Compatible with Mercurius, Carpo and Ceres models. Contact your dealer for up to date technical information.

Impeccable design



The Movo canopy is made of premium quality high strength polycarbonate and cast aluminum to minimize weight and maximize durability.



Its completely molded manufacturing process provides superior design, quality and stability.



The fabrics and zippers used are water repellent, robust and of high quality. The materials have been selected to ensure the durability of the product, while protecting its user.



The windows are made from nautical grade 3 UV resistant mica, which optimizes visibility and the durability of the canopy.

movo
evolution

Distributed by the VERMEIREN Group



Developed in Canada

Developed and designed by Quebec engineers.



Super-fast installation



Ultralight



Superior quality

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **MOV0 Raincover** complies with the highest norms, laid down by certified test institutes in most European countries.



Technical data

■ Total width	736 mm
■ Total height	965 - 1117 mm
■ Total length	1245 mm
■ Folded length	787 mm
■ Interior height	711 mm
■ Interior width	711 mm

■ Interior depth	1346 mm
■ Materials	mica - nautical quality grade 3 anti-UV cast aluminum and high-resistance polycarbonate
■ Washable	yes
■ Weight	8,7 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com

