

2217

FREEDOM

WHITE

The Freedom is a sturdy tricycle ideal for demanding adolescents available in black and white. A perfect tricycle for strolling through countryside.

Range of options available.
Also available as e-Bike.



total length

1599



total width

670



saddle height

720-860



height saddle-pedal

590-780



step-in height

174



height steering

986-1174



wheel size

20" x 1,75



total weight

25,5



users weight

110

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



FREEDOM BLACK

- Different options available.
- The steering wheel (different steerings available) is adjustable in height and angle.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- V brakes at the front with parking lock. Optional disc brake for the rear wheels available.
- Big metal transport tray at the rear.



2217

Inside leg
from 68 to 88 cm

14 years
and
above



V-brakes