

Forest 3

Endurance, Comfort and Adaptability

The **FOREST 3** is the perfect companion for users who are looking for a compact indoor wheelchair and a powerful outdoor model. Endurance, comfort and adaptability are the main characteristics of the Forest 3.



Option SE43
50° adjustment of the
seat and the backrest



Option SE42
Lift function
up to 300 mm



Standard controller



LCD controller



Forest 3

Ideal for indoor and outdoor use. Thanks to the pendulum axle, the Forest 3 is always in contact with the ground, even when the terrain is uneven (up to 120 mm difference).



Option SE55
Reinforced
joystick



Option SE83
Mini joystick



Option B66MAL
Hemiplegic
armrest

Colors (available in matt or glossy)



Technical data

Width	595 / 610 / 650 mm
Length	1185 / 1205 / 1224 mm
Height	1036 mm
Seat height	430 - 480 mm
Seat width	400 / 450 / 500 mm
Seat depth	420 - 520 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Forest 3** complies with the highest norms, laid down by certified test institutes in most European countries.



Option SE39 - Electrical adjustment of the seat and back

Batteries	2 x 12V / 80Ah - 220W or 350W
Speed	6 km/h - 10 km/h - 14 km/h
Driving range	37 km - 50 km
Weight	139 kg
Maximum user weight	130 kg
Additional information	

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeireplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer