

# Quadri-Light

ROLLATOR



total width	650 mm
seat depth	200 mm
seat width	450 mm
seat height	610 mm
total height	805 - 930 mm
total length	700 mm
total weight	7,25 kg
max. users weight	130 kg
wheels (front/rear)	(8") 200 x 40 mm (7") 178 x 40 mm
colours	metallic red - metallic silver

Measurement tolerances:  $\pm 1,5$ mm / 1,5kg.  
We reserve the right to introduce technical changes.

Model shown only for illustration purposes. Certain features may be optional.

**N.V. Vermeiren N.V.** - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: +32 (0)3 620 20 20 • Fax: +32 (0)3 666 48 94 • info@vermeiren.be • website: www.vermeiren.com



Opening up new horizons in mobility