

The logo features a large, stylized number '8' in a gold color with a black outline. Inside the two loops of the '8' are two solid black circles. The word 'CARPO4' is written in a bold, white, sans-serif font across the center of the '8'. Below the '8', the words 'Limited edition' are written in a white, elegant cursive script.

CARPO4  
*Limited edition*



ultimate driving pleasure

# LCD display



new seat  
for a better  
seating comfort



LCD display

exclusive color



new basket



new rims





LCD display

mirror left  
and right

extra strong  
shopping basket

comfortable seat  
with headrest  
and pouch

front and rear lighting  
and indicators

front bumper for  
more security

lateral reflectors  
for a better  
visibility





# technical specifications



- Length : 1450 mm
- Width : 660 mm
- Height : 1200 mm
- Ground clearance : 150 mm
- Batteries : 2x12 V - 75Ah
- Driving range : 45 km\*
- Speed : 15 km/h\*
- Users weight : 150 kg
- Color : brown
- Display : LCD
- Odometer



option: transport box

\* Typical range, calculated under best conditions in accordance with the European Standards. \* Speeds and actionradius can deviate (ca.15%).



NEW Carpo 4 Champagne



NEW Carpo 4 White



Expected: Carpo 4 Deluxe

# NEW ARRIVALS

## N.V. Vermeiren N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium



Tel.: +32 (0) 3 620 20 20  
Fax: +32 (0) 3 666 48 94

For more information, contact your local dealer

[www.vermeiren.com](http://www.vermeiren.com) - [info@vermeiren.be](mailto:info@vermeiren.be)