


690


Tailor made ramps




690 - KIT 1



690 - KIT 2





690 - KIT 3



690 - KIT 4

Our 690 ramp models are the most versatile and sturdy ramps on offer. Adaptable to almost any situation and easy to modify. Determine the height to be traversed, then select the correct ramp kit. With this modular system from 18 mm to 199 mm you can adjust the kit to the desired height and load up to 200 kg.

Available in 5 different sizes from small to extra large.



690 - KIT 5






Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility

Technical data Ramps - 690

■ KIT 1 - length x width x height	250 x 750 x 18-55 mm	■ Weight	1,3 kg
■ KIT 2 - length x width x height	500 x 750 x 36-91 mm	■ Weight	4,3 kg
■ KIT 3 - length x width x height	750 x 750 x 72-127 mm	■ Weight	9,2 kg
■ KIT 4 - length x width x height	1000 x 750 x 108-163 mm	■ Weight	15,8 kg
■ KIT 5 - length x width x height	1250 x 750 x 144-199 mm	■ Weight	20 kg
■ Maximum load	200 kg	■ Additional information	CE

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.




ECO-friendly

All these ramps are made from environment-friendly PE. Non-toxic PE, the same material as found, e.g., on the inside of milk cartons.

PE products are acid-resistant and corrosion-free and resistant to all kinds of household detergents.

Your guarantee of quality

Regarding quality, our ramps comply with the highest norms, laid down by certified test institutes in most European countries.



N.V. VERMEIREN N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility



Subject to technical alterations