


Navix

rearwheeldrive














optional cushions L34, L35 and BZ7-E elektronisch adjustable legrests













Pictures for illustration purposes only. Actual product may differ from the images shown.

Opening up new horizons in mobility

Technical information Navix rearwheeldrive

 total width	 seat width	 seat inclination	 back inclination	 height armrests	 seat height	 total height	 total length	 seat depth	 back height	 max. users weight
580	400	5°/7,5°/10°	10°/20°/30°/40°	260-350	400/450/500	995	1010	380-500	580	130
580	450	5°/7,5°/10°	10°/20°/30°/40°	260-350	400/450/500	995	1010	380-500	580	130
630	500	5°/7,5°/10°	10°/20°/30°/40°	260-350	400/450/500	995	1010	380-500	580	130

Deviations in sizes and weight (15 mm and 1,5 kg) are possible

 weight	 batteries	 power	 speed	 driving range	 max. slope upwards	 max. slope downwards	 ground clearance	 turning radius	 wheels
79	2x12V 38Ah	2 x 220 W	6 km/h	20 km	9°	9°	60	1250	rear: air 317,5x57,15 front: PU 200x50

Class A wheelchair for indoor use.

Options

- Plexi Table
- LED Lights
- Anti-tipping
- Inclinable footplate
- Sideward controller
- Different seat cushions
- Different back cushions
- Side support
- Headrest
- Electrical legrests
- Controller adjustable in height

Colors



LIFT

300 mm lift for a better SOCIAL INTEGRATION

The lift allows me to interact with people around me, enhancing my mobility and social integration.

300 mm lift for a better ACCESSIBILITY

- Kitchen cupboards
- Cash dispensers
- Bar
- etc...



N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN
Opening up new horizons in mobility

Subject to technical alterations