


Albatros

The Vermeiren Albatros was designed to bring a real comfort to the patient during the transfers but also for the people who are taking care of the patient.

The Albatros is intended for people with walking difficulties or no walking abilities.



Completely foldable for transport..



Optional: slings

Control box with a red emergency stop and an internal battery charger



Caster with brake

SLINGS	
Available in M - XL (solid)	
M	XL
40 - 120 kg	100 - 200 kg

Use only Vermeiren approved accessories.



ALBATROS	1040	580	940	1040 - 1730	500	44,8	150
----------	------	-----	-----	-------------	-----	------	-----

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations.



N.V. Vermeiren N.V. - Vermeirenplein 1-15 - 2920 Kalmthout - Belgium
 Tel: 00 32(0)3 620 20 20 • Fax: 00 32(0)3 666 48 94
 www.vermeiren.com • info@vermeiren.be