

# 286i

## ROLLATOR

Width	600 mm
Length (unfolded)	680 mm
Seat width	420 mm
Seat height	560 mm
Brakes	Multifunctional brake
Max. inclination	use on flat underground
Weight	9 kg
With basket & tray	9,85 kg
Maximum user weight	130 kg
Maximum load basket	5 kg
Wheels	Ø 190 mm x 45 mm
Color	blue

Subject to technical alterations. Measurements in mm, weight in kg.  
Deviations in sizes and weight (1,5 mm and 1,5 kg) are possible.



**N.V. Vermeiren N.V.** - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: +32 (0)3 620 20 20 - Fax: +32 (0)3 666 48 94 - info@vermeiren.be - website: www.vermeiren.com  
Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility