

D200P

Standard wheelchair

The wheelchair **D200P** is an excellent combination between technicality, comfort and **VERMEIREN** quality.




Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility

Technical data D200P

Total width	580 / 600 / 620 / 640 / 660 / 680/700/720 mm
Folded width	300 mm
Seat width	380 / 400/ 420 / 440 / 460 / 480/500/520 mm
Height armrests	220 - 240 mm
Seat height	470 - 530 mm
Height backrest	420 mm

Total height	880 - 940 mm
Seat depth	430 mm
Total length	1010 mm
Weight	16 kg
Users weight	130 kg
Additional info	

Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible

Colours



C71
blue



Options

The **D200P** can be equipped as transit model (T30) or extended with several options such as travelwheels, anti-tipping or a B05 armrest. Please consult your dealer or the orderform for a complete and actual list.

The wheelchair must be adapted to you, not vice versa

Our basic principle is: the wheelchair must be adapted to you, not vice versa. This is why our wheelchairs have a modular design.

All types of wheelchairs can be adapted, have numerous accessories and all options are mutually interchangeable. This is the result of our continuous research and more than 64 years of experience.

Your guarantee of quality

Regarding quality, D200P complies with the highest norms, laid down by certified test institutes in most European countries.



B78 - Anti-tipping wheels



B86 - Travelwheels



B05 - Adjustable armrest



Transit model (T30)

N.V. VERMEIREN N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility



Subject to technical alterations