

# 9062 XXL

## XXL TOILET CHAIR

XXL commode with a seat width of 600 mm and max. users weight of 175 kg.



**C11**  
black

max. users weight  
**175 kg**



bucket with handle



780	<b>600</b>	500	540	400	15	175
-----	------------	-----	-----	-----	----	-----

# Interval XXL



Metal side rails version




max. users weight  
**270 kg**

The 3-functions Interval XXL bed has electrical adjustments and is available in 1200 and 1400 mm. Delivered in a white frame with a grid mattress base. The legpart is provided with an electrical adjustable knee-buckle whereas the footpart can be adjusted manually.

Max. users weight: 270 kg.



Delivered on a transportset.  
Installed and ready for usage in a few minutes.



1200x2000	1330 1530	2080	430	830	-	30°	70°	16°	275
-----------	--------------	------	-----	-----	---	-----	-----	-----	-----

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

INTERVAL - 9062

