

Interval








Electrically adjustable home care bed



By using a simple control box, this high/low bed can easily be adjusted to the desired height. From the normal access height, the bed can be lowered to ease the problem of lifting the legs onto the bed. The bed can also be raised to ease the user to the standing position. When nursing, the high position eases the burden on the carer.



Technical information Interval

 width	 length	 min. position	 max. position	 inclination head support	 inclination knee buckle	 weight
1030	2080	360	790	83°	30°	95 kg

Sizes in mm, weight in kg

Deviations in sizes and weight (1,5 mm and 1,5 kg) are possible

Comfortable and practical

The Interval bed has an electrically adjustable head support. By simply pressing a button, the head support moves into a more vertical position, enabling you to read, eat or watch television. The footrest can be placed in the correct position manually.

This high/low bed has 4 wheels, each with a separate brake, therefore guaranteeing stability. The strong and flexible coated baseboard is easy to maintain. The bed can be assembled and disassembled quickly and without tools.

Aesthetic and easy to maintain

The wooden front and rear panels give this bed a homely look. The warm wooden appearance makes the bed an attractive piece of furniture, suitable for all surroundings and easy to maintain.

Options

- Side panels in beech
- Bed mounted self-raising device: removable and can be fitted onto the right or left corner of the bed

Adaptable to your needs

For the past 55 years, Vermeiren has played an important role in home care products. We continually strive to meet your individual needs.

Our products are the result of continuous research and many years of experience.

N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN
Opening up new horizons in mobility