


EP00211 / I

HEAD SIDE WITH WOODEN LATHS				Plan No.
LUNA 2   LUNA BASIC 2 (997) 1/2				EP00211 / I
Scale: 1:6	Drafted	Verified	Confirmed	
Date	08/01/2016	13/07/2022	13/07/2022	
Initials	KR	PM	AZ	

Material			
No.	Plan	Article	Quantity
1	YP05077DP	2899337	1
2	YP05229	2704837	1
3	YP08733	2626300	2
4	X001016	1001682	2
5	M8x40DIN933	2993705	2
6	XP01421	2993116	7
7	XP01420	2993117	14
8	---	1001637	2
9	D8,4 DIN125	1003199	4
10	---	1998885	2
11	---	---	-
12	M8X40	1008739	4
13	---	1902971	1
14	---	1995837	1
15	---	1905446	1
16	Motor LA27	1996201	1
17	---	1996203	1
18	---	1996204	1
19	---	1996205	2
20	---	1005899	2
21	---	1996206	1
22	X004967	1908262	2
23*	---	W620228	1

\* - not depicted / nie ujeta

EP00211 / I

HEAD SIDE WITH WOODEN LATHS			Plan No.
LUNA 2   LUNA BASIC 2 (997) 2/2			EP00211 / I
Scale: 1:6	Drafted	Verified	Confirmed
Date	08/01/2016	13/07/2022	13/07/2022
Initials	KR	PM	AZ

