


2991041


V200 DRIVE SYSTEM HEM2				Plan No. E000651 / G
Sheet 1/2				
Scale: 1:6	Drafted	Verified	Confirmed	
Date	29/09/2010	11/12/2023	11/12/2023	
Initials	DD	AZ	TK	

E000651 / G

Material			
No.	Plan	Article	Quantity
1	X002976	1995343	2
2A	P00065423	2684100	1
2B	P00069096	2684300	1
3	Y017173	1742454 (R)	1
4	Y017205	1742554 (L)	1
5	M6x40 ISO7380	1006111	4
6	M6 DIN986	1008492	4
7	---	1901208	1
8	---	1901745	1
9	---	1994568	4
10	DIN471	1995259	4
11A	P00064457	2683500	1
11B	P00071833	2683600	1
11C	P00071834	2683700	1
12A	X005892	G103100	1
12B	X005891	G103000	1
13	X005893	1915848	1
15	M6 DIN7349	1003321	1
16	M6x25	1903394	1
17	P00064455	GP29900	2
18	D3x26 DIN 1481	1908821	2
19	5x5x10 DIN6885A	2998153	2

2991041

2A Seat width 38 to 52
 2B Seat widths 53 to 60
 11A Seat widths 38 to 44
 11B Seat widths 46 to 52
 11C Seat widths 53 to 60
 12A Seat widths 38 to 44
 12B Seat widths 46 to 60

V200 DRIVE SYSTEM HEM2				Plan No.
Sheet 2/2				E000651 / G
Scale: 1:6	Drafted	Verified	Confirmed	
Date	29/09/2010	11/12/2023	11/12/2023	
Initials	DD	AZ	TK	

E000651 / G