



1910521

| | | | | |
|------------------------|------------|------------|------------|---|
| SEAT FRAME L35 FOREST3 | | | | Plan No. E000953 / J |
| Sheet 1/2 | | | | |
| Scale: 4:25 | Drafted | Verified | Confirmed |  |
| Date | 16/04/2014 | 31/08/2021 | 31/08/2021 | |
| Initials | SD | SD | PJ | |

E000953 / J

| Material | | | |
|----------|---------------------|--------------|----------|
| No. | Plan | Article | Quantity |
| 1 | X001016 (M8 DIN985) | 1001682 | 6 |
| 2 | X003403 | 1617400 | 4 |
| 3 | Y017414 | 1775254 (40) | 1 |
| | Y017984 | 1769754 (45) | 1 |
| | Y017889 | 1776054 (50) | 1 |
| 4 | Y018411 | 1778354 (40) | 1 |
| | Y018368 | 1778454 (45) | 1 |
| | Y018369 | 1778554 (50) | 1 |
| 4* | Y024071 | 2080754 (40) | 1 |
| | Y024008 | 2078754 (45) | 1 |
| | Y024072 | 2080854 (50) | 1 |
| 5 | Y019775 | 1790554 | 1 |
| 6 | M6 DIN985 | 1001680 | 12 |
| 7 | M6x12 | 1007000 | 4 |
| 8 | M12 PA6.6 | 1993362 | 2 |
| 9 | M8x20 DIN7991/8.8 | 1903136 | 4 |
| 10 | M8x40 | 1008739 | 2 |
| 11 | M6x30 DIN7984/8.8 | 1992721 | 12 |
| 13 | 30x15 | 1009992 | 2 |
| 14 | Y021021 | 2007454 | 1 |
| 15 | Ø4x11 | 1001583 | 4 |

| Material | | | Seat width | | | Quantity |
|----------|------------|---------------------------|------------|---------|---------|----------|
| No. | Seat depth | Colour | 40 | 45 | 50 | |
| 12 | 42-48 | Black skai / Blue lining | 1908874 | 1908875 | 1908876 | 1 |
| | 42-48 | Black skai / Green lining | 1908868 | 1908869 | 1908870 | 1 |
| | 42-48 | Black skai / Red lining | 1908871 | 1908872 | 1908873 | 1 |
| | 42-48 | Black skai / Grey soft | 1907467 | 1907468 | 1907469 | 1 |
| | 50-52 | Black skai / Blue lining | 2995195 | 2995196 | 2995197 | 1 |
| | 50-52 | Black skai / Green lining | 2995198 | 2995199 | 2995200 | 1 |
| | 50-52 | Black skai / Red lining | 2995201 | 2995202 | 2995203 | 1 |
| | 50-52 | Black skai / Grey soft | 2995204 | 2995205 | 2995206 | 1 |

* for BZ10 attachment

1910521

E000953 / J

| | | | | |
|------------------------|------------|------------|------------|---|
| SEAT FRAME L35 FOREST3 | | | | Plan No. |
| Sheet 2/2 | | | | E000953 / J |
| Scale: 4:25 | Drafted | Verified | Confirmed |  |
| Date | 16/04/2014 | 31/08/2021 | 31/08/2021 | |
| Initials | SD | SD | PJ | |