



NO.	Material		
	Plan	Article	Number
1	YC00921	1C322	1
	YC00922	1C323	1
2	XP00092	2990042	1
	XP00093	2990043	1
3	X003371	1997774	2
4	X003368	1997773	2
5	-	1007000	4
6	-	1990752	12
7	-	1007954	4
8	DIN 912 M6X10	1001491	2
9	-	1552000	2
10	-	3990082	2

WHEEL CHAIR 101 ARMREST B02

Plan No. / Plannr.

EC00485/A

Scale / Schaal: 1 : 5

Drafted / Getekend

Verified / Nagezien

Confirmed / Bevestigd

Date / Datum

01/06/2020

01/06/2020

01/06/2020

Init. / Init.

AL

JJQ

TG

