



Bill of Material			
No.	Plan	Article	Quantity
1	Side frame	E00050067	1
2	Cross	E00050163	1
3	Rear wheels	E00051594	1
4	Front wheels	E00051691	1
5	Backrest	E00049624	1
6	Sports brake	E00050162	1
	Push/Pull brake	E00050534	1
7	B06	E00048657	1
	B06 IN	EC00219	1
	*B06 / 90°	E00087278	1
	BA6	EC00305	1
	BK6	E00050417	1
	B08	E000981	1
	BK7/BL7	EP00499	1
8	B82(fender)	E00050783	1
	Bi5/Bi5A	E00050770	1
9	Upholstery	E00048650	1
	Upholstery options	E00051931	1
10	B78	E00045457	1
11	B90	E000734	1

\* - use only with front wheel`s size 4" or 5"

V500 ACTIVE EXPLODED			
Scale: 1:8	Drafted	Verified	Confirmed
Date	2023-02-03	22/08/2024	22/08/2024
Initials	AS	AZ	AS

--
Plan No.
E00051729/B1.1